



Your Trust, Our Comfort.

SENNY DENTAL LIMITED

Company Profile

Senmy Dental: Your *Trusted* Dental Chair Partner

✓ Who We Are

We are a leading dental unit manufacturer with over 10 years of export experience.

We offer high-quality dental chairs and a wide range of dental equipments and accessories, providing a one-stop solution for dental clinics and distributors.

✓ Our Advantages

Strict quality control to ensure reliability

Excellent service and fast response

Support for OEM & customized orders

✓ Global Reach

Our products are exported to 50+ countries across Europe, America, the Middle East, and Asia. We've built long-term partnerships with local dealers worldwide.

✓ Let's Connect

Looking for a trusted dental supplier? Contact us now to explore long-term cooperation.

About Us



S106



S103Pro



S101Pro



S101



S102



S102Plus



S102Pro



S103



S103Plus



S104



S104Pro



S104Plus






S105Pro



S105Plus





Dental Chair



S109



S109Pro



S109Plus



S109L

Dentist Stool




S301T



S302T

Color Card

Soft Leather



Micro Fiber
Soft Leather



LED Sensor Lamp

L01



L05



L09B Dental Unit Type



L03



L03Pro



L09C
Ceiling-Rack Type



L06



L07



L08



Dental Suction Unit



Suction Unit



Suction Unit



Movable Suction Unit



Movable Suction Unit

Dental Aircompressor



Air Compressor
for 1-2 Chairs



Air Compressor
for 3-4 Chairs



Air Compressor
for 5-6 Chairs




Air Compressor
for 8 Chairs



Air Compressor
for 10 Chairs

Accessories



Endoscope System
Wifi Mode



Top Mounted



Dental Handpiece Kit



Built-in Light Curing



Built-in LED
Ultrasonic Scaler



Built-in Teeth
Whitening Unit



Built-in Electric Motor



Selection Valve
for Suction






Built-in Fiber
Optic System



3-Way Syringe

Structural Reference



Configuration of Dental Unit



Configuration	S101	S101Pro	S102	S102Plus	S102Pro	S103	S103Pro	S103Plus	S104	S104Pro	S104Plus	S105	S105Pro	S105Plus
24V silent DC motor (TAIWAN TIMOTION)	○	○	○	○	○	○	●	●	○	●	●	○	○	○
Auto spittoon flush and cup filler control system	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Assistant control system	○	○	●	○	●	●	●	●	●	●	●	●	●	●
Rotatable spittoon	—	●	●	●	●	●	●	●	●	●	●	●	●	●
Saliva ejector with strong and weak suction	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Double articulated headrest	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Purified water supply system	●	●	●	●	●	●	●	●	●	●	●	●	●	●
3-way syringe (cold/hot water)	●	●	●	●	●	●	●	●	●	●	●	●	●	●
LED film viewer	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Multi-function foot pedal	○	○	○	○	●	●	●	●	●	●	●	●	●	●
LED operating lamp with sensor	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Dentist stool	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Soft leather cushion	●	●	●	●	●	●	○	○	●	○	○	○	○	○
Micro fiber soft leather cushion	○	○	○	○	○	○	●	●	○	●	●	●	●	●
Low mounted instrument tray	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Top mounted instrument tray	○	○	○	○	○	○	○	○	—	—	—	○	—	—
Ultrasonic scaler	○	○	○	○	○	○	○	○	○	○	○	○	○	○
LED curing light	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Endoscope system	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Electric motor	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Handpiece fiber optic system	○	○	○	○	○	○	○	○	○	○	○	○	○	○

● : Standard Equip

○ : Options

— : No-design



Senmy Dental Limited

Mobile/Wechat/Whatsapp: +86 18937180915

Email: info@senmydental.com

Web: www.senmydental.com

Add: Danzao Town, Nanhai District, Foshan City, Guangdong Province, China.



Web



Whatsapp