



Senmy



Dental Autoclave

Making dentistry safer and better.

—— www.senmydental.com ——





Dental Autoclave

E101(Class B)	8L	12L
Voltage/Power	110-240V/1800W	110-240V/1800W
Chamber inner dimensions (mm)	Φ170×320	Φ200×368
Unit dimensions (mm)	566×408×425	566×408×425
Package size (mm)	700×540×550	700×540×550
Depth (cm)	60	60
Net weight/Grossweight(kg)	40/44	42/47



Dental Autoclave

E101(Class B)	18L	23L
Voltage/Power	110-240V/1800W	110-240V/1800W
Chamber inner dimensions (mm)	Φ250×360	Φ250×450
Unit dimensions (mm)	625×438×473	625×438×473
Package size (mm)	770×580×570	770×580×570
Depth (cm)	60	60
Net weight/Grossweight(kg)	52/55	55/58



Dental Autoclave

E102(Class N)	23L
Voltage/Power	110-240V/1800W
Chamber inner dimensions (mm)	Φ250×450
Unit dimensions (mm)	625×438×473
Package size (mm)	770×580×570
Depth (cm)	60
Net weight/Grossweight(kg)	42/47





E103 (Class N) Auto Drain

Chamber Capacity	18L
Chamber Size	Ø250XD350mm
Autoclave Outer Size	56cm*39cm*36cm
Sterilize Temp	122°C & 132°C
Sterilize Pressure (Mpa)	0.12Mpa & 0.23Mpa
Sterilize Time	20Min & 6Min
Material	Stainless Steel SUS#304



E103Pro (Class N) Auto Water Supply

Chamber Capacity	18L
Chamber Size	Ø250XD350mm
Autoclave Outer Size	56cm*39cm*36cm
Sterilize Temp	122°C & 132°C
Sterilize Pressure (Mpa)	0.12Mpa & 0.23Mpa
Sterilize Time	20Min & 6Min
Material	Stainless Steel SUS#304



E104 UV Sterilizer

Capacity	27L
Voltage	220V
Frequency	50 Hz
Power	15 W
Sterilization Time	15min
Volume	36cm*28cm*59cm
Weight	8.9kg



E105 Portable Pressure Steam Sterilizer

Capacity	18L
Voltage	220V 50 Hz
Rated Working Pressure	0.14~0.16MPa
Max. Safe Pressure	0.165MPa
Power of the Electro-thermal Tube	2kW
Dimension	41cm*41cm*43cm
Weight	16kg



Sterilizer Pouches

5.5cm*200m	12rolls/carton
7.5cm*200m	8rolls/carton
10cm*200m	6rolls/carton
15cm*200m	4rolls/carton
20cm*200m	4rolls/carton
.....	



E301 Dental Sealing Machine

Power	AC 220/110 50Hz
Voltage	500W
Fuse pipe	3TA
Width of Sealing	250mm(max)
Width of pressure bar	12mm
Gross weight	8 KGS
Net weight	5.5 KGS
Size	430*340*250(mm)



E302 Water Distiller

Capacity	4L
Rated Voltage	220V
Rated Frequency	50 Hz
Rated Power	800 W
Distilling Quantity	1L/H
Over-heat Fuse	blows at 160°C
Material	Stainless Steel Inner Tank and Cover



E303 Ultrasonic Cleaner

Cleaning Tank Volume	2.5L
Power Cable Size	≤1.5m
Ultrasonic Power	70W
Ultrasonic Frequency	40KHz
Tank's Dimension	247*150*78mm
Time Setting	60/120/180/240/300/360/480S

**One-press heating up to temperature of 65°C.

E305 Handpiece Lubricator

Volume	300ML
Voltage	110V/220V
Power	40W
Air Pressure	4kg/cm ² -8kg/cm ²
Working Pressure	2.5kg/cm ²
Working Temperature	5-40degree
Net Weight	8.5KG
Carton Size	40*46*48cm



Senmy Dental Limited

Danzao Town, Nanhai District, Foshan City, Guangdong Province, China.
Mobile/Whatsapp/Wechat: +86 18937180915



Scan the QR code
for more information