



LED SENSOR LAMP



Professional Dental Supplier

Mobile: +86 18937180915
Email: info@senmydental.com
Web: www.senmydental.com
Add: Danzao Town, Nanhai District, Foshan City,
Guangdong Province, China



whatsapp



wechat

Senmy Dental Limited

L01

LED口腔灯

LED Dental Lamp

Senmy




Appearance Patent No

外观专利号:

201430388365.7

双色温

Double color temperature



● 技术参数

输入电压/Input Voltage	AC / DC 12-24V
功率/Power	12W
光照度/Illumination	8000-15000Lx
色温/Color Temperature	4500K / 3000K
光斑尺寸/Spot size	80 × 160mm(Distance for 700mm)
LED灯数量/Amount of LED tube	4 × 2个/pcs
包装尺寸/Dimension of carton	39.5 × 34.5 × 20cm
包装重量/净重 GW/NW	2.1KG /1.4KG
每包装数量/Amount of one carton	1台/piece

● Specifications

L02

LED口腔灯


LED DENTAL LAMP



Appearance Patent No

外观专利号:

202230515211.4



进口LED

● 技术参数

输入电压/Input Voltage

功率/Power

光照度/Illumination

色温/Color Temperature

光斑尺寸/Spot size

LED灯数量/Amount of LED tube

包装尺寸/Dimension of carton

包装重量/净重 GW/NW

每包装数量/Amount of one carton

● Specifications

AC / DC 12-24V

7W

8000-25000Lx

4500K

80 × 160mm(Distance for 700mm)

4个/pcs

33 × 33 × 14cm

1.3KG/1KG

1台/piece

L03

LED口腔灯

LED Dental Lamp



Senmy



● 技术参数

输入电压/Input Voltage	AC / DC 12-24V
功率/Power	5W
光照度/Illumination	8000-25000Lx
色温/Color Temperature	5000K Color Temperature
光斑尺寸/Spot size	80 × 160mm (Distance for 700mm)
LED灯数量/Amount of LED tube	1个/pes
包装尺寸/Dimension of carton	43.5 × 36.5 × 19cm
包装重量/净重 GW/NW	1.7KG / 1KG
每包装数量/Amount of one carton	1台/piece

● Specifications

输入电压/Input Voltage	AC / DC 12-24V
功率/Power	5W
光照度/Illumination	8000-25000Lx
色温/Color Temperature	5000K Color Temperature
光斑尺寸/Spot size	80 × 160mm (Distance for 700mm)
LED灯数量/Amount of LED tube	1个/pes
包装尺寸/Dimension of carton	43.5 × 36.5 × 19cm
包装重量/净重 GW/NW	1.7KG / 1KG
每包装数量/Amount of one carton	1台/piece

L04

LED口腔灯

LED Dental Lamp



Senmy



● 技术参数

输入电压/Input Voltage

功率/Power

光照度/Illumination

色温/Color Temperature

光斑尺寸/Spot size

LED灯数量/Amount of LED tube

包装尺寸/Dimension of carton

包装重量/净重 GW/NW

每包装数量/Amount of one carton

● Specifications

AC / DC 12-24V

5W

8000-25000Lx

5000K Color Temperature

80 × 160mm(Distance for 700mm)

1个/pcs

43.5 × 36.5 × 19cm

1.7KG / 1KG

1台/piece

L05

LED口腔灯

LED DENTAL LAMP

Senmy



Appearance Patent No

外观专利号:

202230055547.7

进口LED

色温可调

触摸调节



● 技术参数

输入电压/Input Voltage

功率/Power

光照度/Illumination

色温/Color Temperature

光斑尺寸/Spot size

LED灯数量/Amount of LED tube

包装尺寸/Dimension of carton

包装重量/净重 GW/NW

每包装数量/Amount of one carton

● Specifications

AC / DC 12-24V

9W

8000-30000Lx

3500-5500K 可调Adjust Color Temperature

80 × 160mm(Distance for 700mm)

6个/pcs

38 × 36 × 21.5cm

2.8KG / 1.8KG

1合/piece

L06

LED口腔灯

LED Dental Lamp

Senmy



● 技术参数

输入电压/Input Voltage	AC12V-AV24V
功率/Power	15W
光照度/Illumination	8000-25000Lx
色温/Color Temperature	3000K、4500K、6000K
光斑尺寸/Spot size	80 × 160mm(Distance for 700mm)
LED灯数量/Amount of LED tube	6个/pes
包装尺寸/Dimension of carton	38x36x 21.5cm
包装重量/净重 GW/NW	2.8KG/2.3KG
每包装数量/Amount of one carton	1台/piece

● Specifications

L06Pro

LED口腔灯

LED Dental Lamp



● 技术参数	● Specifications
输入电压/Input Voltage	AC12V-AV24V
功率/Power	15W
光照度/Illumination	8000-25000Lx
色温/Color Temperature	3000K、4500K、6000K
光斑尺寸/Spot size	80 × 160mm(Distance for 700mm)
LED灯数量/Amount of LED tube	8个/pcs
包装尺寸/Dimension of carton	38x36x 21.5cm
包装重量/净重 GW/NW	2.8KG/2.3KG
每包装数量/Amount of one carton	1台/piece

L07

LED口腔灯

LED DENTAL LAMP



- Imported LED
- Adjustable Color Temperature
- New - generation IMD Process Panel

● 技术参数	● Specifications
输入电压/Input Voltage	AC 12-24V
功率/Power	9 VA
光照度/Illumination	8000-30000Lx
色温/Color Temperature	3500-5500K 可调Adjust Color Temperature
光斑尺寸/Spot size	80 × 160mm(Distance for 700mm)
LED灯数量/Amount of LED tube	6个/pcs
包装尺寸/Dimension of carton	43 × 30 × 20cm
包装重量/净重 GW/NW	2.2KG/1.5KG
每包装数量/Amount of one carton	1台/piece



L08

LED口腔灯

LED Dental Lamp



Senmy

- Imported LED
- Adjustable Color Temperature
- New - generation IMD Process Panel



● 技术参数

输入电压/Input Voltage	AC 12-24V
功率/Power	9 VA
光照度/Illumination	8000-30000Lx
色温/Color Temperature	3500-5500K 可调Adjust Color Temperature
光斑尺寸/Spot size	80 × 160mm(Distance for 700mm)
LED灯数量/Amount of LED tube	6个/pcs
包装尺寸/Dimension of carton	42 × 39 × 22cm
包装重量/净重 GW/NW	2.8KG / 1.9KG
每包装数量/Amount of one carton	1台/piece

● Specifications

L09B Built-in Type

LED口腔种植灯 LED DENTAL PLANTING LAMP



Appearance Patent No

外观专利号:

201830109868.4

- Lens Magnification: 2 Million
- Magnification Rate: 20
- Optical Zoom Magnification: 10
- Image Resolution: 1080p
- Monitor: 20" 1080p HDMI
- Storage Capacity: 10GB
- Wireless Network: 2.4G/5.8G Dual Frequency
- Voice Control: Support
- LAN Screen Projection: Support
- Live Streaming: Support
- Remote Video: Support
- Online Updates: Support

Light Source: Imported LED
Lights: 26 LEDs
Lamp Size: $\Phi 380\text{mm}$
Color Temperature: 4000-5500k Adjustable
Illumination: 30,000-90000LX
Activity Adjustment: 360° Adjustable



内置术野高清摄录头
Built-in HD camcorder (Optional)

L09A Moving Type

LED口腔种植灯 LED DENTAL PLANTING LAMP



L09C Ceiling Rack

LED口腔种植灯 LED DENTAL PLANTING LAMP

Senmy



● 技术参数	● Specifications
输入电压/Input Voltage	AC / DC 12-24V
功率/Power	30W
光照度/Illumination	30000-90000Lx
色温/Color Temperature	4000-5500K 可调Adjust Color Temperature
光斑尺寸/Spot size	Φ120mm(Distance for 900mm)
LED灯数量/Amount of LED tube	26个/pcs
包装尺寸/Dimension of carton	60 × 50 × 28cm
包装重量/净重 GW/NW	7KG / 4.6KG
每包装数量/Amount of one carton	1台/piece